

Vascular Accessories



Doppler Rollcart CT-1

For use with ES-100X Surgical MiniDop®, Bidop® 3 and ES-100VX MiniDoppler®

CT-1 \$328.00



Doppler Rollcart CT-2

For use with HD307 and Smartdop® 45

CT-2 \$465.00



Doppler Rollcart CT-30 (for clinic use)

Holds the Doppler, laptop and accessories for convenient vascular testing with Smart-V-Link® software and a compatible Doppler. For use with Smartdop® 30EX, Smartdop® 45, Smartdop® 50EX and Bidop® 3

CT-30 \$295.00



Vascular Cuffs

Reusable vascular cuffs are high quality and exceptionally durable. VC-10 and VC-12 cuffs are latex and rubber free. Velcro fastening and high density cuff stitching ensure a secure fit during inflation.

DVC-1.9 (1.9 cm) digit cuff	\$42.00
DPC-2.5 (2.5 cm) digit/penile cuff	\$42.00
VC-7.5 (7.5 cm) transmetatarsal cuff	\$42.00
VC-10 (10 cm) arm/ankle cuff	\$42.00
VC-12 (12 cm) lower extremity cuff	\$42.00



Footpedal Remote (FP-1)

Convenient footpedal remote control for use with the MiniDop® ES-100X and HD307 surgical Dopplers. Enables user operation without assistance in the sterile field.

FP-1 \$186.00



Speaker (XSP-1)

External speaker for use with the ES-100X MiniDop®. Ideal when additional sound volume is desired. Clips to Doppler rollcart (CT-1) for user convenience.

XSP-1 \$49.00



External Charger with 4 AA Rechargeable Batteries (RAA-1)

External battery charger and 4 AA rechargeable batteries for use with the ES-100X Surgical MiniDop®. Save money on battery replacements.

RAA-1 \$38.00



AC Adaptor

AC adaptor/charger for Smartdop® series and HD307 Dopplers.

Adaptor/Charger \$45.00

Vascular Accessories



Doppler & Accessories Bag (CS-1)

Lightweight padded black case with shoulder strap and handle to easily transport Doppler and accessories. Includes a separate compartment to fit a laptop computer.

CS-1 \$59.00



Headset (HD-1)

Offers user convenience in high ambient noise conditions or when private sound assessment is preferred. The main Doppler speaker will cut off when the headset is in use.

Headset \$22.00



Sterile



Non-Sterile

Ultrasound Gel

Sterile, Surgical (0101) 20 gm foil pack

48 packs per box \$94.00

Non-Sterile (0102) 60 gm tube

12 tubes per box \$28.00



Hand Held Sphygmomanometer (SPG-1)

Aneroid bedside bulb and gauge with button controlled adjustable deflate mechanism. Use with vascular cuffs.

SPG-1 \$85.00



Disinfectant Spray

Use to safely and effectively clean and disinfect Dopplers and probe transducers (clinic only)

Disinfectant Spray 12 oz. \$14.95



Printer Paper

Smardop® series Dopplers (PA-1)

Box of 5 rolls \$38.00

VASCULAR INSTRUMENTS
KOVEN
 TECHNOLOGY INCORPORATED
Koven Technology, Inc.
 Your Vascular Healthcare Partner

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com