

Diagnosis is the First Step

Choose the Bidop® 3 Vascular Doppler for early diagnosis of vascular disease


- Fully customizable - choose **only** the options & accessories to meet the specific needs of your facility
- Real-time waveform LCD - displays the waveform during testing
- 30 waveform memory storage - save waveforms to view later or download to Smart-V-Link® software for EHR's, DICOM® compatible files & documentation for reimbursement
- Automatic gain and baseline control - no manual adjustments required
- Flat vascular probe available - eliminates the need to hold the probe during testing
- PPG available to expand vascular testing capabilities
- Waveform FREQUENCY ANALYSIS available when used with optional Smart-V-Link®
- Doppler rollcart available for user convenience
- Spanish display option



Web based training available for unlimited access to on-demand training & support



*Convenient, handheld design ...
Results you can see*



Smart-V-Link® vascular testing software is available for Waveform Frequency Analysis, EHR's, DICOM® compatible files & hard-copy reports for reimbursement without usage fees or significant hardware costs (run on any Windows® based computer).

GLN: 0894912002005

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com



Find us on Facebook