

# *Sensitive to your Surgical Doppler Demands*

## *Choose the ES-1000SPM Surgical Smartdop™*



- *Wide selection of interchangeable probe styles and frequencies for user flexibility*
- *Autoclavable probes available*
- *Specialty probes for Vascular, Laparoscopic, Plastic & Reconstructive, Cardiovascular, and Neurosurgery*
- *Surgical and non-invasive probes available in 5, 8, 10 and 20 MHz*  
*(non-invasive available in 5, 8 & 10 MHz)*
- *1-year warranty on Doppler and probes reduces replacement costs*  
*(See warranty information for CRP, NRP, and ACP series probes)*
- *Compact and battery operated*
- *Integrated printer for documentation*
- *CSA Approved*

*Lightweight & Portable ...  
Document Results*



**Pencil Probe**  
(P8M05S8A 8 MHz)  
\* also available in 5 & 10 MHz



**Surgical Specialty Probes**  
(8, 10, & 20 MHz available)  
Cardiovascular, Laparoscopic, Neurovascular  
and autoclavable models available



**Non-invasive Probes**  
(5, 8, and 10 MHz available)

VASCULAR INSTRUMENTS  
**KOVEN**  
TECHNOLOGY CANADA INC.

**Koven Technology Canada Inc.**  
*Your Vascular Healthcare Partner*

T: 204.943.2596 / 204.943.2624 • F: 204.943.2652 • [info@koven.com](mailto:info@koven.com) • [www.koven.ca](http://www.koven.ca) • email PO's: [canadapurchase@koven.com](mailto:canadapurchase@koven.com)