


Diagnosis is the First Step

Choose the Smartdop™ 45 Vascular Testing System for early diagnosis of vascular disease

- **Convenient hand-held size includes features found in larger systems**
- **Fully customizable - choose only the options & accessories to meet the specific needs of your facility**
- **Real-time waveform LCD - displays the waveform before printing**
- **30 waveform memory storage - save waveforms to view later or download to Smart-V-Link® software for EHR's & DICOM® compatible files**
- **Automatic gain and baseline control - no manual adjustments are required**
- **Waveform strip generated for immediate documentation**
- **Flat vascular probe available - Eliminates the need to hold the probe during testing**
- **PPG and PV modules available to expand testing capabilities**
- **French display option**



Web based training available for unlimited access to on-demand training & support



Smart-V-Link® vascular testing software is available for EHR's & DICOM® compatible files without usage fees or significant hardware costs (run on any Windows® based computer).

*Unique,
handheld
design...
Proven results*

VASCULAR INSTRUMENTS
KOVEN
 TECHNOLOGY CANADA INC.
 Koven Technology Canada Inc.
 Your Vascular Healthcare Partner

T: 204.943.2596 / 204.943.2624 • F: 204.943.2652 • info@koven.com • www.koven.ca



Find us on Facebook