

Sensitive to your Surgical Probe Demands

- **Durable and Cost-effective**
- **Surgical probes are interchangeable with Koven Surgical Dopplers**
- **Specialty probes for Vascular, Cardiovascular, Laparoscopic, Neurosurgery, Plastic and Reconstructive**
- **Pre-sterile single use & autoclavable probes available**
- **Wide variety of sterilization options available (ETO, Steris®, Steris® V-Pro, Sterrad®, Optreoz® 125L, Autoclave)**
- **Non-invasive probes also available (2, 5, 8, and 10 MHz)**
- **1 Year warranty ***

* 1 year parts & labor warranty when used as recommended. Excludes VDP-8 single-use. Warranties for ACP, NRP, and CRP series probes differ. See warranty information for details.



Surgical Pencil Probes

Available in 5, 8, & 10 MHz

Sterilization: ETO, Sterrad, Optreoz 125L



Autoclavable Probe

Model: ACP-8 (8 MHz)

Sterilization: Autoclave

* Adaptor required



Surgical Thin Probes

Models: VRP-8 & VRP-20 (8 & 20 MHz)

Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro, Optreoz® 125L

* Adaptor required



Single-Use, Pre-Sterile Probes

Model: VDP-8 (8 MHz)

* Adaptor required



Cardiovascular Probes

Models: CRP-10H & CRP-20H (10 & 20 MHz)

Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro, Optreoz® 125L

* Adaptor required



Bayonet Style Neurovascular Probe

Models: NRP-10H & NRP-20H (10 & 20 MHz)

Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro, Optreoz® 125L

* Adaptor required



Laparoscopic Probe

Model: LRP-8 (8 MHz)

Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro, Optreoz® 125L

* Adaptor required



Flat Monitoring Probe

Model: FDP-8 (8 MHz)

Sterilization: ETO, Sterrad®, Steris® V-Pro

* Adaptor required



Air Emboli Flat Monitoring Probe

Model: BF2M20S8A (2 MHz)

Compatible with ES-100VX and Bidop® 3 Noninvasive Dopplers

Sterilization: Sterrad®, Optreoz® 125L

