

Hands-Free Vascular Testing

by Koven

Choose the *Smartdop® 30EX Vascular Testing System* for hands-free testing & diagnosis of vascular disease


Leasing options available



*Durable & Effortless...
Automated, reliable
vascular testing*



Web based training available for unlimited access to on-demand training & support



Smart-V-Link® vascular testing software is available for EHR's & DICOM® compatible files without usage fees or significant hardware costs (run on any Windows® based computer).

- **Use with the flat vascular probe or PPG probe to eliminate the need to hold the probe during testing** (pencil style probes also available)
- **Makes vascular testing with Doppler faster & easier - automatically calculates ABI and TBI ratios**
- **Fully customizable - choose only the options & accessories to meet the specific needs of your facility**
- **Includes features previously found in larger, more costly systems**
- **Built-in automatic cuff inflator quickly performs :**
ABIs (10 min.)
Segmental Pressures (15 min.)
PVRs (5 min.)
PPG Toe Pressures
- **30 waveform memory storage - save patient study results to view later or download to Smart-V-Link® software for EHR's & DICOM® compatible files**
- **Reimbursable vascular studies add revenue to your practice**
- **Doppler rollcart available for user convenience**
- **Spanish display option**

GLN: 0894912002005

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner



Find us on Facebook

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com